

Product Datasheet

FOXL2 Antibody (orb1671903)

Description	FOXL2 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, EM, IHC, IP, WB
Immunogen	Amyloid-beta peptide.
Target	FOXL2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P05067
Expiration Date	12 months from date of receipt.