


Product Datasheet

ITGB7 Antibody (orb1671938)

Description	ITGB7 Antibody
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC, IP
Immunogen	FIB27 was prepared by immunizing rats with TK1 cells.
Target	ITGB7
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P26011, P26010
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed hum...