


Product Datasheet

CD19 Antibody (orb1671950)

Description	CD19 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Human prolymphocytic leukaemia cell line JVM3.
Target	CD19
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2a kappa
Clonality	Recombinant
Uniprot ID	P15391
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat sple...