



Product Datasheet

Prion Antibody [3F4], Mouse IgG1 (orb1671972)

Description	Prion Antibody [3F4], Mouse IgG1
Reactivity	Hamster, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Hamster 263K PrPs.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed U25...



Western Blot using anti-major Prion prot...