


Product Datasheet

CD55 Antibody (orb1671977)

Description	CD55 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Human CD55.
Target	CD55
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG lambda
Clonality	Recombinant
Uniprot ID	P08174
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lu...



Flow-cytometry using the anti-CD55 (Comp...