



Product Datasheet

C3 Antibody (orb1671987)

Description	C3 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC, WB
Immunogen	Human EGFR.
Target	C3
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P00533
Expiration Date	12 months from date of receipt.



Flow-cytometry
using the anti-EGFR
S3 12...



Immunofluorescence
staining of fixed
NIH...