

Product Datasheet

NTRK1 Antibody (orb1671998)

Description	NTRK1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, IHC
Immunogen	Human TrkA ECD.
Target	NTRK1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P04629
Expiration Date	12 months from date of receipt.