


Product Datasheet

GCGR Antibody (orb1672014)

Description	GCGR Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, IF
Immunogen	Recombinant human glucagon receptor.
Target	GCGR
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P47871
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed Hep...