


Product Datasheet

PXN Antibody (orb1672020)

Description	PXN Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide corresponding to paxillin-LD2 motif (residues 141-156) fused to MBP; peptide sequence: NLSELDRLLELNAVQHNP.
Target	PXN
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P49023
Expiration Date	12 months from date of receipt.



Western Blot
using anti-
Paxillin (LD2
Mo...