



Product Datasheet

Cd274 Antibody (orb1672031)

Description	Cd274 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Murine PDL1.
Target	Cd274
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q9EP73
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-PDL1 antibody ...



Flow-cytometry using anti-PD-L1 antibody...