


Product Datasheet

SLC3A2 Antibody (orb1672034)

Description	SLC3A2 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Epitope peptide (LLHGDFHAFSAG-GGGCCD) linked to KLH.
Target	SLC3A2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P08195, Q01650
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed HeL...