



Product Datasheet

SNCG Antibody (orb1672053)

Description	SNCG Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The KLH-conjugated peptide "YLAPSGPSGTLK", derived from Gamma-synuclein. Gamma-synuclein plays a role in neurofilament network integrity.
Target	SNCG
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q6FHG5
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-SNCG antib...



Immunofluorescence staining of fixed U25...