



Product Datasheet

AMHR2 Antibody (orb1672061)

Description	AMHR2 Antibody
Reactivity	Human
Conjugation	Unconjugated
Immunogen	The KLH-conjugated peptide "SFVLNLGK" derived from Galectin 1. Galectin 1 binds beta-galactoside and other complex carbohydrates and may regulate apoptosis, cell proliferation and cell differentiation.
Target	AMHR2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P09382
Expiration Date	12 months from date of receipt.



Western Blot using anti-LGALS1 antibody ...



Flow cytometry using the Anti-LGALS1 ant...