



Product Datasheet

BCL5 Antibody (orb1672066)

Description	BCL5 Antibody
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Peptide "DGVVEITGK" derived from Heat shock protein beta-1, involved in stress resistance and actin organisation, conjugated to KLH.
Target	BCL5
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P04792
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...



Flow cytometry using the Anti-HSPB1 anti...



Western Blot using anti-HSPB1 antibody S...