



Product Datasheet

SEPTIN5 Antibody (orb1672078)

Description	SEPTIN5 Antibody
Reactivity	Human, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Immunoprecipitate from human brain.
Target	SEPTIN5
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	Q99719
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rin usin...



Western Blot using anti-CDCrel-1 antibod...