


Product Datasheet

FCGR2B Antibody (orb1672091)

Description	FCGR2B Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human CD32.
Target	FCGR2B
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P31994
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-CD32 antibody ...