

Product Datasheet

CD7 Antibody (orb1672109)

Description	CD7 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC
Target	CD7
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2b kappa
Clonality	Recombinant
Uniprot ID	P09564
Expiration Date	12 months from date of receipt.