



Product Datasheet

Polysialic acid Antibody [735], Rabbit IgG (orb1672113)


Description	Polysialic acid Antibody [735], Rabbit IgG
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP
Immunogen	Homopolymers of α 2,8-linked sialic acid of Escherichia coli K1.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of guinea p...



Immunohistochemical staining of rat spin...



Western Blot using anti-polysialic acid ...