



Product Datasheet

CD4 Antibody (orb1672132)


Description	CD4 Antibody
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, IP, WB
Target	CD4
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P01730
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-CD4 antibody M...



Flow-cytometry using anti-CD4 antibody M...



Flow-cytometry using anti-CD4 antibody M...