


Product Datasheet

NOTCH2 Antibody (orb1672169)

Description	NOTCH2 Antibody
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, IF
Immunogen	NRR of mouse Notch 2.
Target	NOTCH2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	O35516, Q04721
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed HeLa...