


## Product Datasheet

### Siglec1 Antibody (orb1672192)

|                            |                                                                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Siglec1 Antibody                                                                                                                                                                       |
| <b>Reactivity</b>          | Mouse                                                                                                                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                           |
| <b>Tested Applications</b> | IHC, WB                                                                                                                                                                                |
| <b>Immunogen</b>           | Spleen cells from AO rats immunised with thioglycollate-elicited murine peritoneal macrophages (TPM) that had been cultured in mouse serum to induce sheep erythrocyte receptor (SER). |
| <b>Target</b>              | Siglec1                                                                                                                                                                                |
| <b>Preservatives</b>       | PBS with 0.02% Proclin 300.                                                                                                                                                            |
| <b>Concentration</b>       | batch dependent                                                                                                                                                                        |
| <b>Storage</b>             | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.                                                                                                       |
| <b>Note</b>                | For research use only                                                                                                                                                                  |
| <b>Isotype</b>             | IgG2a kappa                                                                                                                                                                            |
| <b>Clonality</b>           | Recombinant                                                                                                                                                                            |
| <b>Uniprot ID</b>          | <a href="#">Q62230</a>                                                                                                                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                        |



Anti-CD169  
staining of  
murine  
macrophage...