


Product Datasheet

TNF Biosimilar Antibody (orb1672199)

Description	TNF Biosimilar Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, In vivo
Immunogen	Purified recombinant human TNF (rTNF).
Target	TNF
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P01375
Expiration Date	12 months from date of receipt.



Western Blot using anti-TNF alpha antibo...