


Product Datasheet

ABCB1 Antibody (orb1672200)

Description	ABCB1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IP
Immunogen	Drug resistant K-562 cells.
Target	ABCB1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P08183
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-P-glycopro...