


Product Datasheet

FGA Antibody (orb1672247)

Description	FGA Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Human cross-linked D-Dimer.
Target	FGA
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P02671
Expiration Date	12 months from date of receipt.



Binding curve of anti-D-Dimer antibody 3...