



Product Datasheet

CD40 Antibody (orb1672250)


Description	CD40 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Activation, FC, IHC, IP
Immunogen	Human erythrocyte-rosette-negative tonsillar lymphocytes.
Target	CD40
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P25942
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-CD40 antibody ...



Flow-cytometry using anti-CD40 antibody ...



Flow-cytometry using anti-CD40 antibody ...