

Product Datasheet

CD38 Antibody (orb1672253)

Description	CD38 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Viable cells (1 x 10E7) of the human cell line, BJAB (Human Burkitt lymphoma B cell line), which was were injected intraperitoneally into 8-week-old BALB/c mice.
Target	CD38
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P28907
Expiration Date	12 months from date of receipt.