



Product Datasheet

CD52 Antibody (orb1672258)

Description	CD52 Antibody
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	Cytotoxic Assay, ELISA, FC, IHC-Fr, IHC-P, In vivo, WB
Immunogen	Human lymphocytes.
Target	CD52
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P31358
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-CD52 antibody ...



Flow-cytometry using anti-CD52 antibody ...