



Product Datasheet

CD52 Antibody (orb1672259)


Description	CD52 Antibody
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	Cytotoxic Assay, ELISA, FC, IHC-Fr, IHC-P, In vivo, WB
Immunogen	Human lymphocytes.
Target	CD52
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P31358
Expiration Date	12 months from date of receipt.



Immunofluorescent staining of rat thymus...



Immunohistochemical staining of human th...



Immunofluorescence staining of fixed rat...