



Product Datasheet

FABP3 antibody (orb167226)


Description	Rabbit polyclonal antibody to FABP3.
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human FABP3(NP_004093.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14KDa
Uniprot ID	P05413
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...



Immunohistochemical staining of human st...



Immunofluorescence analysis of U2OS cell...