


Product Datasheet

CD3E Antibody (orb1672262)

Description	CD3E Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human peripheral E+ cells.
Target	CD3E
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P07766
Expiration Date	12 months from date of receipt.



Immunofluorescence
staining of fixed
Mol...