

Product Datasheet

Fluorescein Antibody [4-4-20 (enhanced)], Human IgG2 (orb1672308)

Description	Fluorescein Antibody [4-4-20 (enhanced)], Human IgG2
Conjugation	Unconjugated
Tested Applications	Control
Immunogen	Fluorescein I-bovine serum albumin (FII-BSA).
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2 kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.