

Product Datasheet

His tag (C-term) Antibody [3D5], Mouse IgG2b (orb1672319)

Description	His tag (C-term) Antibody [3D5], Mouse IgG2b
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified his tagged Fab fragment.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2b kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.