



Product Datasheet

KLKB1 antibody (orb167232)


Description	Rabbit polyclonal antibody to KLKB1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human KLKB1(NP_000883.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80KDa
Uniprot ID	P03952
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of rat kidn...



Immunohistochemical staining of rat brai...