


## Product Datasheet

### His tag (C-term) Antibody [3D5], Mouse IgG1 (orb1672327)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | His tag (C-term) Antibody [3D5], Mouse IgG1                                      |
| <b>Conjugation</b>         | Unconjugated                                                                     |
| <b>Tested Applications</b> | ELISA, WB                                                                        |
| <b>Immunogen</b>           | Purified his tagged Fab fragment.                                                |
| <b>Preservatives</b>       | PBS with 0.02% Proclin 300.                                                      |
| <b>Concentration</b>       | batch dependent                                                                  |
| <b>Storage</b>             | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| <b>Note</b>                | For research use only                                                            |
| <b>Isotype</b>             | IgG1 kappa                                                                       |
| <b>Clonality</b>           | Recombinant                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                  |



orb1672327  
(1 ug/ml -  
right panel)  
stain...