



Product Datasheet

MYC Antibody (orb1672344)

Description	MYC Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, IP, WB
Immunogen	A synthetic peptide corresponding to aa 408-439 from C-terminus of human c-myc.
Target	MYC
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P01106
Expiration Date	12 months from date of receipt.



Left Panel:
orb1672344
(2 ug/ml)
stainin...



ELISA of anti-
epitope tag
antibodies on

...