



Product Datasheet

IFI16 Antibody (orb1672345)


Description	IFI16 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Anti-IFI16 antibody (orb1672345) was raised against a peptide corresponding to 16 amino acids near the carboxyl terminus of human IFI16. The immunogen is located within the last 50 amino acids of IFI16.
Target	IFI16
Preservatives	IFI16 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	IFI16 antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 88kD Observed: 88 kD
Uniprot ID	Q16666
NCBI	Q16666.3
Expiration Date	12 months from date of receipt.



WB
Validation in
Human and
Rat Cell
Line...



Western Blot
Validation in
Mouse
Tissues...



Western Blot
Validation in
Rat Testis.
L...