

## Product Datasheet

### Tonsil Membrane Lysate (orb1672362)

|                            |                                                                                                                                                                                                                                      |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Tonsil Membrane Lysate                                                                                                                                                                                                               |
| <b>Tested Applications</b> | WB                                                                                                                                                                                                                                   |
| <b>Preservatives</b>       | HEPES pH 7.9, MgCl <sub>2</sub> , KCl, EDTA, Sucrose, Glycerol, Sodium deoxycholate, and NP-40.                                                                                                                                      |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                                                                               |
| <b>Concentration</b>       | 2 mg/mL or better                                                                                                                                                                                                                    |
| <b>Storage</b>             | Store at 2-8°C for continuous use. For extended storage, freeze working aliquots at -70°C. Repeated freezing and thawing is not recommended. Under proper storage conditions the shelf life is half a year from the date of receipt. |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                |
| <b>Application notes</b>   | Tonsil membrane protein lysate is for use in Western blotting, 10 µg to 20 µg per lane is recommended for mini gel.                                                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                      |