

## Product Datasheet

### Stomach-Cardia Membrane Lysate (orb1672374)

<b>Description</b>	Stomach-Cardia Membrane Lysate
<b>Tested Applications</b>	WB
<b>Preservatives</b>	HEPES pH 7.9, MgCl <sub>2</sub> , KCl, EDTA, Sucrose, Glycerol, Sodium deoxycholate, and NP-40.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	2 mg/mL or better
<b>Storage</b>	Store at 2-8°C for continuous use. For extended storage, freeze working aliquots at -70°C. Repeated freezing and thawing is not recommended. Under proper storage conditions the shelf life is half a year from the date of receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	Stomach: Cardia membrane protein lysate is for use in Western blotting, 10 µg to 20 µg per lane is recommended for mini gel.
<b>Expiration Date</b>	12 months from date of receipt.