



Product Datasheet

RAD17 antibody (orb167273)


Description	Rabbit polyclonal antibody to RAD17
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human RAD17(NP_002864.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	77kDa
Uniprot ID	075943
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of MCF-7 cel...



Immunohistochemical staining of human CO...