


Product Datasheet

CHMP2B antibody (orb167313)

Description	Rabbit polyclonal antibody to CHMP2B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human CHMP2B (NP_054762.2).
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	28kDa
Uniprot ID	Q9UQN3
Entrez	25978
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...