

Product Datasheet

Recombinant Human MEC (CCL28) (orb1673593)

Description	Recombinant Human MEC (CCL28)
Form/Appearance	Recombinant CCL-28 was lyophilized from PBS pH 7.4.
Concentration	1.0 EU/ μ g of recombinant protein as determined by the LAL method.
Storage	Upon reconstitution, this cytokine can be stored in working aliquots at 2° - 8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95% as determined by SDS-PAGE
Source	E. coli
MW	12 kDa
Solubility (25°C)	A quick spin of the vial followed by reconstitution in distilled water to a concentration not less than 0.1 mg/mL. This solution can then be diluted into other buffers.
Expiration Date	6 months from date of receipt.