



Product Datasheet

MAGEA1 antibody (orb167376)


Description	Rabbit polyclonal antibody to MAGEA1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human MAGE-1/MAGE1A (NP_004979.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	45kDa
Uniprot ID	P43355
Entrez	4100
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...



Immunohistochemical staining of human br...



Immunohistochemical staining of human li...