




Product Datasheet

TXNRD2 antibody (orb167396)

Description	Rabbit polyclonal antibody to TXNRD2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human Thioredoxin reductase 2 (Thioredoxin reductase 2 (TXNRD2)) (NP_006431.2).
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	56kDa
Uniprot ID	Q9NNW7
Entrez	10587
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...

Immunohistochemistry
of paraffin-embedde...Immunohistochemistry
of paraffin-embedde...