



Product Datasheet

UMPS antibody (orb167398)


Description	Rabbit polyclonal antibody to UMPS.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human UMPS(NP_000364.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56kDa
Uniprot ID	P11172
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human colon tissue.



Immunohistochemical staining of human colon tissue.



Immunofluorescence analysis of HeLa cells.