



Product Datasheet

SNAIL antibody (orb167414)


Description	Rabbit polyclonal antibody to SNAIL
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human Snail(NP_005976.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29kDa
Uniprot ID	O95863
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human co...



Immunohistochemical staining of mouse br...



Immunohistochemical staining of rat brai...