

## Product Datasheet

### Recombinant Mouse Dexamethasone-induced Ras-related protein 1 (Rasd1) (orb1674328)

**Description**

Recombinant Mouse Dexamethasone-induced Ras-related protein

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Endotoxins</b>        | Not test.                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Preservatives</b>     | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.                                                                                                                                                                          |
| <b>Form/Appearance</b>   | Liquid or Lyophilized powder                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                             |
| <b>Tag</b>               | N-terminal 10xHis-tagged                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Application notes</b> | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Protein Sequence</b>  | MKLAAMIKKMCPSDSELSIPAKNCYRMVILGSSKVGKTAIVSRFLTGRFEDAYTPTI<br>EDFHRKFYSIRGEVYQLDILDTSGNHPPFAMRRLSILTGDVFILVFSLDNRDSFEE<br>VQRLKQQLDTSCLKNKTKENVDPVPLVICGNKGDRDFYREVEQREIEQLVGDD<br>PQRCAYFEISAKKNSSLDQMFRALFAMAKLPSEMSPDLHRKVSQYCDVLHKKA<br>LRNKKLLRAGSGGGGDHGDAGILAPFARRPSVHSDLMYIREKTSVGSQAKDKE<br>RCVIS                                                                                                   |
| <b>Purity</b>            | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>MW</b>                | 34.2 kDa                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Uniprot ID</b>        | <b>035626</b>                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                                                                                                                                                                                                                                                              |