

Product Datasheet

**Recombinant Severe acute respiratory syndrome
coronavirus 2 Replicase polyprotein 1ab (rep)
(orb1674352)**

Description

Recombinant Severe acute respiratory syndrome coronavirus 2

Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Protein Sequence	AVGACVLCNSQTSRLRCGACIRRPFLCCKCCYDHVISTSHKLVLSVNPYVCNAPGC DVTDTVTLQLYLGGMSSYYCKSHKPPISFPLCANGQVFGLYKNTCVGSDNVTDFNAI ATCDWTNAGDYILANTCTERLKLFAAETLKATEETFKLSYGIATVREVLSRELHL SWEVGKPRPPLNRNYVFTGYRVTKNSKVQIGEYTFEKG DYGD VVYRGTTTYKL NVGDYFVLTSHTVMPLSAPTLVPQEHYVRITGLYPTLNISDEFSSNVANYQKVG QKYSTLQGPPTGKSHFAIGLALYPSARIVYTACSHAAVDALCEKALKYLPIDKC SRIIPARARVECFDKFVNSTLEQYVFCTVNALPETTADIVVFDEISMATNYDLSVV NARLRKHYVYIGDPAQLPAPRTLLTKGTLEPEYFNSVCRLMKTIGPDMFLGTCR RCPAEIVDTVSAVLDNKLKAHKDKSAQCFKMFYKGVITHDVSSAINRPQIGVVR EFLTRNPAWRKAVFISPYNSQNAVASKILGLPTQTVDSQSEYDYVIFTQTETETA HSCNVNRFNVAITRAKVGILCIMSDDRLYDKLQFTSLEIPRRNVATLQ
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Yeast
MW	92.3 kDa
Uniprot ID	P0DTD1
Expiration Date	6 months from date of receipt.