



Product Datasheet

CFI antibody (orb167455)

Description	Rabbit polyclonal antibody to CFI.
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CFI(NP_000195.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60kDa
Uniprot ID	P05156
Dilution Range	WB: 1:500-2000, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of human placenta



Western blot analysis of extracts of Human serum