

## Product Datasheet

**Recombinant Severe acute respiratory syndrome  
coronavirus 2 Replicase polyprotein 1ab (rep)  
(orb1674575)**

## Description

Recombinant Severe acute respiratory syndrome coronavirus 2

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Endotoxins</b>        | Not test.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Preservatives</b>     | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Form/Appearance</b>   | Liquid or Lyophilized powder                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Tag</b>               | C-terminal 6xHis-tagged                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Application notes</b> | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                                                                               |
| <b>Protein Sequence</b>  | AVGACVLCNSQTSRLRCGACIRRPFLCCKCCYDHVISTSHKLVLNVNPNVCNAPGC<br>DVTDTVTLQLYLGGMSSYYCKSHKPPISFPLCANGQVFGLYKNTCVGSDNVTDFNAI<br>ATCDWTNAGDYILANTCTERLKLFAAETLKATEETFKLSYGIATVREVLSRELHL<br>SWEVGKPRPPLNRNYVFTGYRVTKNSKVQIGEYTFEKG DYGDVAVYRGTTTYKL<br>NVGDYFVLTSHTVMPLSAPTLVPQEHYVRITGLYPTLNISDEFSSNVANYQKVG<br>QKYSTLQGPPTGKSHFAIGLALYPSARIVYTACSHAAVDALCEKALKYLPIDKC<br>SRIIPARARVECFDKFVNSTLEQYVFCTVNALPETTADIVVFDEISMATNYDLSVV<br>NARLRKHYVYIGDPAQLPAPRTLLTKGTLEPEYFNSVCRLMKTIGPDMFLGTCR<br>RCPAEIVDTVSAVLDNKLKAHKDKSAQCFKMFYKGVITHDVSSAINRPQIGVVR<br>EFLTRNPAWRKAVFISPYNSQNAVASKILGLPTQTVDSQSEYDYVIFTQTETA<br>HSCNVNRFNVAITRAKVGILCIMSDDRDLQFTSLEIPRRNVATLQ |
| <b>Purity</b>            | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>MW</b>                | 70.0 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Uniprot ID</b>        | <b>P0DTD1</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |