



Product Datasheet

FKBP4 antibody (orb167475)


Description	Rabbit polyclonal antibody to FKBP4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein of human FKBP52/FKBP52/FKBP4(NP_002005.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52KDa
Uniprot ID	Q02790
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human hy...



Immunohistochemical staining of human ad...



Western blot analysis of extracts of var...